

### CLIMSEAT, COOLING SEAT CUSHION

People are not physiologically programmed to sit for long periods of time. This is why ClimSeat cooling seat cushion will become your best ally in your daily life.



ClimSeat will correct the harmful effects of sitting, such as postural and circulatory disorders, hyperthermia (increased temperature) and musculoskeletal problems due to the pressure of sitting on the spine and coccyx.

### PRINCIPLE

ClimSeat is a cooling seat cushion that will improve your daily comfort, promote endogenous testosterone production, increase male fertility and encourage good back support.

It will also relieve post-operative pain, pain from recent childbirth, circulatory problems (including cellulite) and hemorrhoids.



This is the first product to combine:

- **An orthopedic foam and U-shaped cushion**
- **Innovative interchangeable cooling pockets.**



One bag weighs 1.8 kg and can be used without cooling: it already absorbs excess body heat naturally and effectively. For even more cooling, it can be placed in the fridge.

→ With two cooling bags, you will always have a cool one ready to replace the other one saturated with your body heat.

### USE

ClimSeat is very easy to use.

It is a U-shaped seat cushion with a net on top to hold one of the two cooling pockets. This net has Velcro fasteners to hold the cooling bag in place.

To prevent the seat cushion from slipping from the seat, it has a non-slip surface on the bottom.

The seat cushion should be positioned so that **your tailbone is level with the empty part of the "U"** as shown in the photo.

**If you want the cushion to be cooler**, just place one of the cooling bags (green or blue) in the fridge for about 1 hour. Change the bags as soon as you wish.



## CLEANING AND STORAGE

The seat cushion can be removed with a zip.

If necessary, the cover can be washed in the washing machine at 30°C on a gentle cycle.

The inner foam is not washable. However, you can ventilate it.

The cooling pockets can be washed by hand with a damp sponge over the surface.

Never expose the cooling bags to direct sunlight.

## BENEFITS

- ✓ Mechanical and thermal comfort for daily sitting
- ✓ Lumbar support and spinal decompression
- ✓ Postural improvement
- ✓ Stimulation of testosterone and improvement of male fertility
- ✓ Relief of circulatory disorders and anti-cellulite effect
- ✓ Relief of post-operative and post-delivery situations
- ✓ Relief of hemorrhoids (avoids pressure and promotes circulation)
- ✓ Comfort for pregnant women

### Warranty

*This product is neither a toy nor a medical device.*

*Neo Factory - Climsom cannot guarantee the consequences of the use of ClimSeat.*

*Neo Factory - Climsom guarantees the original purchaser of the product against any manufacturing defect in the product for 2 years from the date of purchase.*

*The guarantee will only result in the replacement of the defective product to the exclusion of any other repair.*



Discover all our products at [www.climsom.com](http://www.climsom.com)

NEO FACTORY / CLIMSOM - 6 allée du commandant Charcot 44000 Nantes - France

+33 (0)2 85 52 44 74 [contact@climsom.com](mailto:contact@climsom.com)